


Anti-Magohb Polyclonal Antibody For Wb, Elisa - Q96A72

Item Number: CM0022907



Introduction

Rabbit polyclonal antibody targeting human MAGOHB for Western blot and ELISA. Validated in human and rat, this antibody supports spliceosome and nonsense-mediated decay research.

[Learn More](#)

Parameter	Value
Product Name	MAGOHB Rabbit pAb
Remarks/Alias	MGN2; mago; magoh; MAGOHB
Species	Human
GeneID (Human)	55110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148 of human MAGOHB (NP_060518.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96A72
Protein Weight	17kDa
Shipping	Ice bag