

Anti-Magp-2/Mfap5 Rabbit Monoclonal Antibody For Wb, Elisa - Q13361

Item Number: CM0005664

Introduction

High-quality rabbit monoclonal antibody targeting human MAGP-2/MFAP5, validated for WB and ELISA. Ideal for studying extracellular matrix biology, hematopoiesis, and cardiovascular research. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	MAGP-2/MFAP5 Rabbit mAb
Remarks/Alias	AAT9; MP25; MAGP2; MAGP-2; MFAP-5; MAGP-2/MFAP5
Species	Human
GeneID (Human)	8076
Immunogen	Recombinant protein of human MAGP-2/MFAP5.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13361
Protein Weight	20kDa
Shipping	Ice bag