

Anti-Magt1 Polyclonal Antibody For Wb, Elisa - Q9H0U3

Item Number: CM0028427



Introduction

Rabbit polyclonal antibody against human MAGT1 (IAP, MagT1), validated for WB and ELISA. Accessory OST complex subunit involved in N-glycosylation and magnesium transport. Ideal for research in glycosylation disorders and immunodeficiencies.

[Learn More](#)

Parameter	Value
Product Name	MAGT1 Rabbit pAb
Remarks/Alias	IAP, XMEN, MRX95, OST3B, CDG1CC, PRO0756, SLC58A1, bA217H1.1, MAGT1
Species	Human
GeneID (Human)	84061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-220 of human MAGT1 (NP_115497.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H0U3
Protein Weight	38kDa
Shipping	Ice bag