

Anti-Major Vault Protein (Lrp/Mvp) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14764

Item Number: CM0006139



Introduction

Rabbit monoclonal antibody targeting human major vault protein (MVP/LRP), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for multidrug resistance and cancer research.

[Learn More](#)

Parameter	Value
Product Name	LRP/MVP Rabbit mAb
Remarks/Alias	LRP; VAULT1; LRP/MVP
Species	Human
GeneID (human)	9961
GeneID	9961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-370 of Human LRP/MVP (NP_034146.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14764
Protein Weight	99kDa
Shipping	Ice bag