

Anti-Major Vault Protein (Mvp) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14764

Item Number: CM0020251



Introduction

Rabbit polyclonal antibody targeting human Major vault protein (MVP/LRP), validated for Western blot, IHC-P, and ELISA applications. Reacts with human and mouse species; suitable for multidrug resistance and nucleocytoplasmic transport studies.

[Learn More](#)

Parameter	Value
Product Name	MVP Rabbit pAb
Remarks/Alias	LRP, VAULT1, MVP
Species	Human
Gene ID (Human)	9961
Gene ID	9961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MVP (NP_005106.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14764
Protein Weight	99kDa
Shipping	Ice bag