

Anti-Mak Polyclonal Antibody For Wb, Ihc-P, Elisa - P20794

Item Number: CM0019913

Introduction

High-quality rabbit polyclonal antibody targeting human MAK (RP62). Validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for ciliary biology, spermatogenesis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAK Rabbit pAb
Remarks/Alias	RP62; MAK
Species	Human
GeneID (Human)	4117
GeneID	4117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-435 of human MAK (NP_005897.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20794
Protein Weight	71kDa
Shipping	Ice bag