


Anti-Mak16 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bxy0

Item Number: CM0021414



Introduction

Rabbit polyclonal antibody targeting human MAK16 homolog (RBM13, NNP78), validated for Western blot and ELISA, with cross-reactivity to mouse and rat. The 35 kDa nucleolar protein suggests involvement in ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	MAK16 Rabbit pAb
Remarks/Alias	RBM13; MAK16L; MAK16
Species	Human
GeneID (Human)	84549
GeneID	84549
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-120 of human MAK16 (NP_115898.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXY0
Protein Weight	35kDa
Shipping	Ice bag