


Anti-Malt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Udy8

Item Number: CM0023844



Introduction

Rabbit polyclonal antibody against human MALT1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for studying NF-kappa-B signaling and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	MALT1 Rabbit pAb
Remarks/Alias	MLT; MLT1; IMD12; PCASP1; MALT1
Species	Human
GeneID (Human)	10892
GeneID	10892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-450 of human MALT1 (NP_776216.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UDY8
Protein Weight	92kDa
Shipping	Ice bag