

Anti-Maml1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92585

Item Number: CM0023369



Introduction

Rabbit polyclonal antibody raised against human MAML1, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Mastermind-like protein 1 is a transcriptional coactivator in Notch signaling.

[Learn More](#)

Parameter	Value
Product Name	MAML1 Rabbit pAb
Remarks/Alias	Mam1; Mam-1; MAML1
Species	Human
GeneID (Human)	9794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MAML1 (NP_055572.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92585
Protein Weight	108kDa
Shipping	Ice bag