

Anti-Maml2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8IzI2

Item Number: CM0024261

Introduction

Rabbit polyclonal antibody against human MAML2, a Notch transcriptional coactivator. Applicable in WB, IF/ICC, and ELISA. Recognizes human, mouse, rat. Suitable for studying Notch signaling and gene regulation.

[Learn More](#)

Parameter	Details
Product Name	MAML2 Rabbit pAb
Remarks/Alias	MAM2; MAM3; MAM-3; MLL-MAML2; MAML2
Species	Human
GeneID	84441 (Human)
Immunogen	Synthetic peptide within amino acids 750-850 of human MAML2 (NP_115803.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZL2
Protein Weight	125kDa
Shipping	Ice bag