

Anti-Man2A1 Polyclonal Antibody For Wb, Elisa - Q16706

Item Number: CM0014908

Introduction

Rabbit polyclonal antibody targeting MAN2A1, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying N-glycan biosynthesis and Golgi function.

[Learn More](#)

Parameter	Value
Product Name	MAN2A1 Rabbit pAb
Remarks/Alias	MANA2; MANII; GOLIM7; AMan II; MAN2A1
Species	Human
GeneID (Human)	4124
GeneID	4124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 865-1144 of human MAN2A1 (NP_002363.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16706
Protein Weight	131kDa
Shipping	Ice bag