

Anti-Man2A2 Polyclonal Antibody For Wb, Elisa - P49641

Item Number: CM0019610

Introduction

Rabbit polyclonal antibody targeting human alpha-mannosidase 2x (MAN2A2). Validated for Western blot and ELISA applications. Suitable for investigating N-glycan processing and Golgi-associated glycosylation pathways.

[Learn More](#)

Parameter	Value
Product Name	MAN2A2 Rabbit pAb
Remarks/Alias	MANA2X; MAN2A2
Species	Human
GeneID (Human)	4122
GeneID	4122
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human MAN2A2 (NP_006113.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49641
Protein Weight	131kDa
Shipping	Ice bag