

Anti-Man2B1 Rabbit Polyclonal Antibody For Wb, Elisa - O00754

Item Number: CM0028034

Introduction

Rabbit polyclonal antibody against human MAN2B1, validated for Western blot and ELISA. Cross-reacts with mouse. Suitable for studying lysosomal alpha-mannosidase function in glycoprotein degradation.

[Learn More](#)

Parameter	Value
Product Name	MAN2B1 Rabbit pAb
Remarks/Alias	MANB; LAMAN; MAN2B1
Species	Human
GeneID (Human)	4125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-850 of human MAN2B1 (NP_000519.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O00754
Protein Weight	114kDa
Shipping	Ice bag