

Anti-Manba Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O00462

Item Number: CM0021430

Introduction

MANBA Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, ELISA in human and mouse samples. Recombinant immunogen. Ideal for studying beta-mannosidase in lysosomal research.

[Learn More](#)

Parameter	Value
Product Name	MANBA Rabbit pAb
Remarks/Alias	MANB1; MANBA
Species	Human
GeneID (Human)	4126
GeneID	4126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-300 of human MANBA (NP_005899.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00462
Protein Weight	101kDa
Shipping	Ice bag