


Anti-Mannose Receptor/Cd206 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P22897

Item Number: CM0008934



Introduction

Recombinant rabbit monoclonal antibody targeting human Mannose Receptor/CD206 (UniProt P22897). Suitable for WB, IF-P, IHC-P, ELISA. Ideal for studying endocytosis, pathogen recognition, and macrophage biology.

[Learn More](#)

Parameter	Value
Product Name	Mannose Receptor/CD206 Rabbit mAb
Remarks/Alias	MMR; hMR; CD206; MRC1L1; CLEC13D; CLEC13DL; bA541I19.1
Species	Human
GeneID (human)	4360
GeneID	4360
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-1389 of human Mannose Receptor/CD206 (NP_002429.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P22897
Protein Weight	166kDa
Shipping	Ice bag