


Anti-Mannose Receptor/Cd206 Polyclonal Antibody For Wb - P22897

Item Number: CM0026864



Introduction

High-quality rabbit polyclonal antibody targeting human Mannose Receptor/CD206 (MRC1), validated for WB, IF-P, and ELISA. Suitable for studies of macrophage endocytosis, pathogen recognition, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	Mannose Receptor/CD206 Rabbit pAb
Remarks/Alias	MMR; hMR; CD206; MRC1L1; CLEC13D; CLEC13DL; bA541I19.1; Mannose Receptor/CD206
Species	Human
Gene ID (Human)	4360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 610-860 of human Mannose Receptor/CD206 (NP_002429.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22897
Protein Weight	166 kDa
Shipping	Ice bag