


Anti-Maob Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27338

Item Number: CM0024100



Introduction

Monoamine Oxidase B (MAOB) Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for neurotransmitter metabolism and neurotoxin studies.

[Learn More](#)

Parameter	Value
Product Name	Monoamine Oxidase B (MAOB) Rabbit pAb
Remarks/Alias	MAOB; Monoamine Oxidase B (MAOB)
Species	Human
GeneID (Human)	4129
GeneID	4129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-490 of human MAOB (NP_000889.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27338
Protein Weight	59kDa
Shipping	Ice bag