

Anti-Map1Lc3A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H492

Item Number: CM0016331



Introduction

High-quality rabbit polyclonal antibody targeting human MAP1LC3A, a key autophagy marker. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for autophagy and cellular stress studies.

[Learn More](#)

Parameter	Value
Product Name	MAP1LC3A Rabbit pAb
Remarks/Alias	LC3; LC3A; ATG8E; MAP1ALC3; MAP1BLC3; MAP1LC3A
Species	Human
Gene ID (Human)	84557
Gene ID	84557
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAP1LC3A (NP_115903.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H492
Protein Weight	14kDa
Shipping	Ice bag