

Anti-Map1Lc3A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H492

Item Number: CM0005137



Introduction

CamelBio provides a rabbit monoclonal anti-MAP1LC3A antibody for autophagy studies, validated in WB, IHC-P, and ELISA. Reacts with human, mouse, rat; detects 14kDa protein. Catalog product for research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	MAP1LC3A Rabbit mAb
Remarks/Alias	LC3; LC3A; ATG8E; MAP1ALC3; MAP1BLC3; MAP1LC3A
Species	Human
Gene ID (Human)	84557
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 42-121 of human LC3A (Q9H492).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H492
Protein Weight	14kDa
Shipping	Ice bag