

Anti-Map1Lc3A Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9H492

Item Number: CM0002243



Introduction

Anti-MAP1LC3A rabbit monoclonal antibody validated for Western blot, immunoprecipitation, and ELISA. Reacts with human, mouse, and rat. Ideal for autophagy research, targeting the autophagosome marker LC3A.

[Learn More](#)

Parameter	Value
Product Name	MAP1LC3A Rabbit mAb
Remarks/Alias	LC3; LC3A; ATG8E; MAP1ALC3; MAP1BLC3; MAP1LC3A
Species	Human
GeneID (Human)	84557
GeneID	84557
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAP1LC3A (Q9H492).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H492
Protein Weight	14kDa
Shipping	Ice bag