

Anti-Map1Lc3C Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bxw4

Item Number: CM0028759



Introduction

Anti-MAP1LC3C rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat MAP1LC3C, a key autophagy-related protein involved in xenophagy and aggrephagy.

[Learn More](#)

Parameter	Value
Product Name	MAP1LC3C Rabbit pAb
Remarks/Alias	LC3C; ATG8J; MAP1LC3C
Species	Human
Gene ID (Human)	440738
Gene ID	440738
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human MAP1LC3C (NP_001004343.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXW4
Protein Weight	17kDa
Shipping	Ice bag