

Anti-Map2 Polyclonal Antibody For Wb, If, Ihc, Elisa - P11137

Item Number: CM0025492

Introduction

Rabbit polyclonal antibody targeting MAP2, suitable for WB, IHC-P, IF-P, and ELISA. Detects human, mouse, and rat MAP2. Ideal for neuroscience studies involving microtubule dynamics and dendritic architecture.

[Learn More](#)

Parameter	Value
Product Name	MAP2 Rabbit pAb
Remarks/Alias	MAP-2; MAP2A; MAP2B; MAP2C; MAP2
Species	Human
Gene ID	4133 (Human)
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1510-1650 of human MAP2 (NP_002365.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11137
Protein Weight	200 kDa
Shipping	Ice bag