

Anti-Map2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11137

Item Number: CM0015426

Introduction

High-quality rabbit polyclonal antibody targeting human MAP2. Recombinant immunogen designed for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. For microtubule and neurobiology studies.

[Learn More](#)

Parameter	Value
Product Name	MAP2 Rabbit pAb
Remarks/Alias	MAP-2; MAP2A; MAP2B; MAP2C; MAP2
Species	Human
Gene ID (Human)	4133
Gene ID	4133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human MAP2 (NP_002365.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11137
Protein Weight	200kDa
Shipping	Ice bag