


Anti-Map2K4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P45985

Item Number: CM0012937



Introduction

Rabbit polyclonal antibody against MAP2K4, validated for knockout (KO) studies. Suitable for WB, IHC-P, ELISA; reacts with human, mouse, rat. Targets the stress-activated MAPK/JNK signaling pathway.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MAP2K4 Rabbit pAb
Remarks/Alias	JNKK; MEK4; MKK4; SEK1; SKK1; JNKK1; SERK1; MAPKK4; PRKMK4; SAPKK1; SAPKK-1; K4
Species	Human
Gene ID (Human)	6416
Gene ID	6416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-399 of human MAP2K4 (NP_003001.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45985
Protein Weight	44kDa
Shipping	Ice bag