


Anti-Map3K1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13233

Item Number: CM0023279



Introduction

Anti-MAP3K1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat. Targets MAP3K1/MEKK1 in ERK/JNK and NF-kB signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MAP3K1 Rabbit pAb
Remarks/Alias	MEKK; MEKK1; SRXY6; MEKK 1; MAPKKK1; K1
Species	Human
GeneID (human)	4214
GeneID	4214
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13233
Protein Weight	164kDa
Shipping	Ice bag