


## Anti-Map3K1 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - Q13233

Item Number: CM0021718



### Introduction

High-quality rabbit polyclonal antibody targeting MAP3K1 (MEKK1) for Western blot, immunofluorescence, immunohistochemistry, and ELISA applications. Cross-reacts with human, mouse, and rat MAP3K1. Suitable for kinase signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                          |
|------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MAP3K1 Rabbit pAb                                                                                                              |
| Remarks/Alias    | MEKK; MEKK1; SRXY6; MEKK 1; MAPKKK1; MAP3K1                                                                                    |
| Species          | Human                                                                                                                          |
| GeneID (Human)   | 4214                                                                                                                           |
| GeneID           | 4214                                                                                                                           |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1400-1512 of human MAP3K1 (NP_005912.1). |
| Source           | Rabbit                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                          |
| Application      | WB, IF/ICC, IF-P, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                              |
| SWISS            | Q13233                                                                                                                         |
| Protein Weight   | 164kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                        |