

Anti-Map3K11 Rabbit Polyclonal Antibody For Wb, Elisa - Q16584

Item Number: CM0020517



Introduction

Rabbit polyclonal antibody targeting human MAP3K11 (MLK3). Suitable for Western blot and ELISA detection. Recognizes mouse MAP3K11. Recombinant protein immunogen (aa 140-260). Useful for studying MAPK pathways and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	MAP3K11 Rabbit pAb
Remarks/Alias	MLK3; PTK1; SPRK; MLK-3; MEKK11; MAP3K11
Species	Human
GeneID (Human)	4296
GeneID	4296
Immunogen	Recombinant Protein (Recombinant fusion protein containing a sequence corresponding to amino acids 140-260 of human MAP3K11, NP_002410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q16584
Protein Weight	93kDa
Shipping	Ice bag