

Anti-Map3K4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6R4

Item Number: CM0019696

Introduction

High-quality rabbit polyclonal MAP3K4 antibody for WB, IHC-P, ELISA applications. Reacts with human and rat MAP3K4 (predicted ~182kDa). Essential tool for studying MAP kinase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	MAP3K4 Rabbit pAb
Remarks/Alias	MTK1; MEKK4; MEKK 4; MAPKKK4; PRO0412; MAP3K4
Species	Human
GeneID	4216
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAP3K4 (NP_005913.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y6R4
Protein Weight	182kDa
Shipping	Ice bag