

Anti-Map3K4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6R4

Item Number: CM0003745

Introduction

This rabbit monoclonal antibody targets human MAP3K4 (MEKK4), a kinase in the MAPK pathway. Validated for WB, IHC-P, and ELISA, it cross-reacts with mouse and rat. Molecular weight 182kDa.

[Learn More](#)

Parameter	Value
Product Name	MAP3K4 Rabbit mAb
Remarks/Alias	MTK1, MEKK4, MEKK 4, MAPKKK4, PRO0412, MAP3K4
Species	Human
GeneID (Human)	4216
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MAP3K4 (Q9Y6R4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6R4
Protein Weight	182kDa
Shipping	Ice bag