

Anti-Map3K5 Polyclonal Antibody For Wb, Elisa - Q99683

Item Number: CM0016450



Introduction

Rabbit polyclonal antibody against MAP3K5/ASK1 (Q99683), validated for WB and ELISA. Detects Human, Mouse, Rat MAP3K5, a key kinase in stress and apoptosis signaling. Ideal for cell signaling and apoptosis research.

[Learn More](#)

Parameter	Specification
Product Name	MAP3K5 Rabbit pAb
Remarks/Alias	ASK1; MEKK5; MAPKKK5; MAP3K5
Species	Human
Gene ID (Human)	4217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human MAP3K5 (NP_005914.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99683
Protein Weight	155kDa
Shipping	Ice bag