

Anti-Map3K8 Rabbit Polyclonal Antibody For Wb, Elisa - P41279

Item Number: CM0015330



Introduction

Rabbit polyclonal antibody targeting human MAP3K8 (COT/TPL2), validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying TLR4-mediated MAPK/ERK signaling, cytokine production, and immune responses.

[Learn More](#)

Parameter	Value
Product Name	MAP3K8 Rabbit pAb
Remarks/Alias	COT, EST, ESTF, TPL2, AURA2, MEKK8, Tpl-2, c-COT, MAP3K8
Species	Human
GeneID (Human)	1326
GeneID	1326
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAP3K8 (NP_005195.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41279
Protein Weight	53kDa
Shipping	Ice bag