

Anti-Map4K2 Rabbit Polyclonal Antibody For Wb & Elisa - Q12851

Item Number: CM0020144

Introduction

MAP4K2 rabbit polyclonal antibody validated for WB and ELISA, recognizing human, mouse, rat. Targets the GCK kinase involved in MAPK and Hippo signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	MAP4K2 Rabbit pAb
Remarks/Alias	GCK; BL44; RAB8IP; MAP4K2
Species	Human
GeneID (Human)	5871
GeneID	5871
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12851
Protein Weight	92kDa
Shipping	Ice bag