

Anti-Map4K3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ivh8

Item Number: CM0027863



Introduction

High-quality rabbit polyclonal antibody targeting human MAP4K3 (GLK/MEKKK3), validated for WB and ELISA. Detects a ~101 kDa band; cross-reacts with mouse and rat. Ideal for Hippo and stress signaling studies.

[Learn More](#)

Parameter	Value
Product Name	MAP4K3 Rabbit pAb
Remarks/Alias	GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1; MAP4K3
Species	Human
GeneID (Human)	8491
GeneID	8491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 755-894 of human MAP4K3 (NP_003609.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IVH8
Protein Weight	101kDa
Shipping	Ice bag