


Anti-Mapk11 (P38B/Sapk2B) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15759

Item Number: CM0028483



Introduction

Rabbit polyclonal antibody against MAPK11 (p38 β /SAPK2b), validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for MAP kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	MAPK11 Rabbit pAb
Remarks/Alias	P38B; SAPK2; p38-2; PRKM11; SAPK2B; p38Beta; P38BETA2; MAPK11
Species	Human
GeneID (Human)	5600
GeneID	5600
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 264-364 of human MAPK11 (NP_002742.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15759
Protein Weight	41kDa
Shipping	Ice bag