

Anti-Mapk13 Polyclonal Antibody For Wb, Elisa - O15264

Item Number: CM0019739



Introduction

Rabbit polyclonal antibody targeting human MAPK13 (p38delta, SAPK4). Validated for Western blot (WB) and ELISA, with cross-reactivity in human, mouse, and rat tissues. Ideal for studying stress-activated kinase signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	MAPK13 Rabbit pAb
Remarks/Alias	SAPK4; PRKM13; MAPK 13; MAPK-13; p38delta; MAPK13
Species	Human
Gene ID	5603
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human MAPK13 (O15264).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15264
Protein Weight	42kDa
Shipping	Ice bag