


Anti-Mapk13 Rabbit Polyclonal Antibody For Wb And Elisa - O15264

Item Number: CM0029292



Introduction

Rabbit polyclonal antibody targeting human MAPK13 (p38δ), suitable for WB and ELISA. Validated for human, mouse, and rat samples. Ideal for studying MAP kinase signaling, stress responses, and cellular differentiation.

[Learn More](#)

Parameter	Value
Product Name	MAPK13 Rabbit pAb
Remarks/Alias	SAPK4; PRKM13; MAPK 13; MAPK-13; p38delta; MAPK13
Species	Human
Gene ID	5603
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-365 of human MAPK13 (NP_002745.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15264
Protein Weight	42kDa
Shipping	Ice bag