

Anti-Mapk1/Mapk3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28482 / P27361

Item Number: CM0021330



Introduction

High-quality rabbit polyclonal antibody against MAPK1 (ERK2) and MAPK3 (ERK1), suitable for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for MAPK/ERK signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	MAPK1/MAPK3 Rabbit pAb
Species	Human
GeneID (human)	5594/5595
GeneID	5594/5595
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human MAPK1/MAPK3 (NP_620407.1/NP_002737.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	36kDa/38kDa/40kDa/41kDa/43kDa
Shipping	Ice bag