


Anti-Mapk1/Mapk3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P28482 / P27361

Item Number: CM0021811



Introduction

Rabbit polyclonal antibody targeting human MAPK1/MAPK3 (ERK1/2). Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for signal transduction studies and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	MAPK1/MAPK3 Rabbit pAb
Remarks/Alias	MAPK1/MAPK3
Species	Human
GeneID (Human)	5594/5595
GeneID	5594/5595
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human MAPK1/MAPK3 (ERK1/2) (P27361/P28482)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	36kDa/41kDa/38kDa/40kDa/43kDa
Shipping	Ice bag