


Anti-Mapk6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16659

Item Number: CM0015239



Introduction

Rabbit polyclonal antibody against human MAPK6 (ERK3), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets atypical MAP kinase involved in MAP2 phosphorylation and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	MAPK6 Rabbit pAb
Remarks/Alias	ERK3; PRKM6; p97MAPK; HsT17250; MAPK6
Species	Human
Gene ID (Human)	5597
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 372-721 of human MAPK6 (NP_002739.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16659
Protein Weight	83kDa
Shipping	Ice bag