

Anti-Mapk8 Rabbit Polyclonal Antibody For Wb, Elisa - P45983

Item Number: CM0022094



Introduction

High-quality rabbit polyclonal antibody against human MAPK8/JNK1, validated in WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studying stress-activated kinase signaling and AP-1 regulation.

[Learn More](#)

Parameter	Value
Product Name	MAPK8 Rabbit pAb
Remarks/Alias	JNK; JNK1; PRKM8; SAPK1; JNK-46; JNK1A2; SAPK1c; JNK21B1/2
Species	Human
GeneID (human)	5599
GeneID	5599
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAPK8 (NP_620637.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P45983
Protein Weight	48kDa
Shipping	Ice bag