

Anti-Mapkap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bpz7

Item Number: CM0026301



Introduction

Rabbit polyclonal antibody against MAPKAP1 (SIN1), recognizing human, mouse, and rat mTORC2 subunit. Validated for WB, IHC-P, ELISA, detects ~59kDa band. Ideal for mTOR signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MAPKAP1 Rabbit pAb
Remarks/Alias	MIP1; SIN1; JC310; SIN1b; SIN1g; MAPKAP1
Species	Human
Gene ID (Human)	79109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-522 of MAPKAP1 (NP_001006618.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPZ7
Protein Weight	59kDa
Shipping	Ice bag