

Anti-Mapkapk-2/Mk2 Monoclonal Antibody For Wb, Elisa - P49137

Item Number: CM0007317



Introduction

Rabbit monoclonal antibody targeting human MAPKAPK-2/MK2. Suitable for Western blot and ELISA; detects human, mouse, and rat MAPKAPK-2/MK2. Ideal for stress signaling, cytokine regulation, and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	MAPKAPK-2/MK2 Rabbit mAb
Remarks/Alias	MK2; MK-2; MAPKAP-K2; MAPKAPK-2/MK2
Species	Human
GeneID	9261
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human MAPKAPK-2/MK2 (NP_116584.2).
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49137
Protein Weight	46kDa
Shipping	Ice bag