

Anti-Mapkapk-2/Mk2 Rabbit Polyclonal Antibody For Wb And Elisa - P49137

Item Number: CM0014470



Introduction

Rabbit polyclonal antibody targeting human MAPKAPK-2/MK2 (P49137). Validated for WB and ELISA with cross-reactivity to rat. Ideal for stress signaling, cytokine regulation, and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	MAPKAPK-2/MK2 Rabbit pAb
Remarks/Alias	MK2; MK-2; MAPKAP-K2; MAPKAPK-2/MK2
Species	Human
GeneID (Human)	9261
GeneID	9261
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human MAPKAPK-2/MK2 (NP_116584.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P49137
Protein Weight	46kDa
Shipping	Ice bag