

Anti-Mapkbp1 Rabbit Polyclonal Antibody For Wb And Elisa - O60336

Item Number: CM0019865

Introduction

Validated for Western blotting and ELISA, this antibody detects MAPKBP1 (JNKBP1) across human, mouse, and rat samples, supporting studies on NOD2 and JNK signaling.

[Learn More](#)

Parameter	Value
Product Name	MAPKBP1 Rabbit pAb
Remarks/Alias	JNKBP1; NPHP20; JNKBP-1; MAPKBP1
Species	Human
GeneID	23005
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAPKBP1 (NP_001122080.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60336
Protein Weight	164kDa
Shipping	Ice bag