

Anti-Mapt Polyclonal Antibody For Wb, If/Icc, Elisa - P10637

Item Number: CM0024668



Introduction

Rabbit polyclonal antibody targeting MAPT/Tau (SWISS: P10637), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat tau isoforms. Suitable for neurodegeneration and microtubule biology studies.

[Learn More](#)

Parameter	Specification
Product Name	MAPT Rabbit pAb
Remarks/Alias	Tau; Mtapt; PHF-tau; MAPT
Species	Mouse
Gene ID (Human)	4137
Gene ID	17762
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3-214 of mouse MAPT (NP_034968.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10637
Protein Weight	35-38 kDa/44 kDa/76 kDa
Shipping	Ice bag