


Anti-March8 Polyclonal Antibody For Wb And Elisa - Q5T0T0

Item Number: CM0014418



Introduction

Rabbit polyclonal antibody against human MARCH8 (E3 ubiquitin-protein ligase MARCHF8), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for innate and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	MARCH8 Rabbit pAb
Remarks/Alias	MIR; CMIR; c-MIR; MARCH8; RNF178; MARCH-VIII
Species	Human
Gene ID (Human)	220972
Gene ID	220972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 232-291 of human MARCH8 (NP_659458.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5T0T0
Protein Weight	33kDa
Shipping	Ice bag