

Anti-March9 Polyclonal Antibody For Wb, Elisa - Q86Yj5

Item Number: CM0014075

Introduction

High-quality rabbit polyclonal antibody targeting human MARCH9 (MARCF9, RNF179, MARCH-IX), an E3 ubiquitin ligase that regulates endocytosis and lysosomal degradation of immune receptors. Suitable for WB and ELISA.

[Learn More](#)

Parameter	Value
Product Name	MARCH9 Rabbit pAb
Remarks/Alias	MARCH9; RNF179; MARCH-IX
Species	Human
Gene ID (Human)	92979
Gene ID	92979
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 277-346 of human MARCH9 (NP_612405.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86YJ5
Protein Weight	38kDa
Shipping	Ice bag