

Anti-Marcks Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P29966

Item Number: CM0015739

Introduction

Rabbit polyclonal antibody targeting human MARCKS, applicable for western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studying actin cytoskeleton dynamics, cell motility, and signal transduction pathways.

[Learn More](#)

Parameter	Value
Product Name	MARCKS Rabbit pAb
Remarks/Alias	MACS; 80K-L; PKCSL; PRKCSL; MARCKS
Species	Human
GeneID (Human)	4082
GeneID	4082
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MARCKS (NP_002347.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29966
Protein Weight	32kDa
Shipping	Ice bag