

Anti-Marcksl1 Rabbit Polyclonal Antibody For Wb And Elisa - P49006

Item Number: CM0023724

Introduction

Rabbit polyclonal antibody against human MARCKSL1 (MARCKS-like protein 1), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studies of actin cytoskeleton and cell migration. 19 kDa.

[Learn More](#)

Parameter	Value
Product Name	MARCKSL1 Rabbit pAb
Remarks/Alias	MARCKSL1; F52; MACMARCKS; MLP; MLP1; MRP; MARCKS like 1
Species	Human
Gene ID	65108
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 95-195 of human MARCKSL1 (NP_075385.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49006
Protein Weight	19kDa
Shipping	Ice bag