

Anti-Marco Polyclonal Antibody For Wb, If-P, Elisa - Q9Uew3

Item Number: CM0015077

Introduction

Rabbit polyclonal antibody against human MARCO (scavenger receptor A6/SCARA2) for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for innate immunity and macrophage research.

[Learn More](#)

Parameter	Value
Product Name	MARCO Rabbit pAb
Remarks/Alias	SR-A6; SCARA2; MARCO
Species	Human
Gene ID (Human)	8685
Gene ID	8685
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 421-520 of human MARCO (NP_006761.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UEW3
Protein Weight	53kDa
Shipping	Ice bag